



SpotCheck Plus detects the presence of glucose and lactose on a surface as a rapid indicator of surface hygiene. Simply swab a surface, release the reagent, and if any residue containing glucose and/or lactose is present, the reagent will turn green. The

more contamination present, the quicker the color change and darker the color. SpotCheck Plus instantly validates the hygiene of a surface, allowing immediate corrective action to be taken if necessary. If it turns green, it isn't clean!


Benefits:

- All-in-one device
- Easy-to-interpret color results
- No instrumentation necessary
- Results in 60 seconds
- Pre-moistened swab for reliable collection, recovery, and detection
- 100% recyclable
- 18 month shelf life at refrigerated temperatures (2-8°C)
- 4 week shelf life at room temperatures (21-25°C)

Sensitivity:


- Detects down to 2.5 µmoles D-glucose or 5.0 µmoles lactose in 60 seconds
- Disregard any color change after 2 minutes

Color Results:



Catalog No.	Description	Quantity
GL100	SpotCheck Plus	100

SpotCheck Plus Procedure:



Related Products:



AllerSnap

AllerSnap is a colorimetric allergen prevention test for broad range allergen screening in 15-30 minutes. AllerSnap detects down to 3 µg protein.



ProClean

ProClean is a quick and easy-to-use device for monitoring the cleanliness of surfaces and safeguarding product quality. ProClean detects protein residues left on a surface after cleaning, giving results in 1-10 minutes without any required instrumentation.



EnSURE Touch

EnSURE™ Touch is an advanced monitoring system that collects, analyzes and reports data from a variety of long-trusted quality tests, providing rapid and accurate sanitation verification data. EnSURE™ Touch features a 5-inch touch screen, wireless sync technology, and cloud-based software.

Catalog No.	Description	Quantity
ALS-100	AllerSnap Allergen Prevention	100
PRO100	ProClean Protein Residue Test	100
ETOUCH	EnSURE Touch	1